



Recombinant Escherichia coli Bifunctional protein GlmU (glmU)

Product Code	CSB-EP484932ENM-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7L878
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain 55989 / EAEC)
Purity	>85% (SDS-PAGE)
Sequence	MLNNAMSVVI LAAGKGTRMY SDLPKVLHTL AGKAMVQHVI DAANELGAAH VHLVYGHGGD LLKQALKDDN LNWVLQAEQL GTGHAMQQAA PFFADDEDIL MLYGDVPLIS VETLQRLRDA KPQGGIGLLT VKLDDPTGYG RITRENGKVT GIVEHKDATD ELRQIQEINT GILIANGADM KRWLAKLTNN NAQGEYYITD IALAYQEGR EIVAVHPQRL SEVEGVNNRL QLSRLERVYQ SEQAEKLLLA GVMLRDPARF DLRGTLTHGR DVEIDTNVII EGNVTLGHRV KIGTGCVIKN SVIGDDCEIS PYTVVEDANL AAACITGPFA RLRPGAELLE GAHVGNFVEM KKARLGKGSK AGHLTYLGDA EIGDNNVIGA GTITCNYDGA NKFKTIIGDD VFVGSQTQLV APVTVGKGAT IAAGTTVTRN VGENALAIR VPQTQKEGWR RPVKKK
Source	E.coli
Target Names	glmU
Protein Names	Recommended name: Bifunctional protein GlmU Including the following 2 domains: UDP-N-acetylglucosamine pyrophosphorylase EC= 2.7.7.23 Alternative name(s): N-acetylglucosamine-1-phosphate uridylyltransferase Glucosamine-1-phosphate
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.