



Recombinant Escherichia coli 8-amino-7-oxononanoate synthase (bioF)

Product Code	CSB-MP484959ENM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LC58
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain 55989 / EAEC)
Purity	>85% (SDS-PAGE)
Sequence	MSWQEKINAA LDARRAADAL RRRYPVAQGA GRWLVADDRQ YLNFSSNDYL GLSHHPQIIR AWQQGAEQFG VGSGGSGHVS GYSVAHQALE EELAEWLGYS RALLFISGFA ANQAVIAAMM AKEDRIVADR LSHASLLEAA SLSPSQLRRF VHNDVTHLAR LLASPCPGQQ LVVTEGVFSM DGDIAPLAEI QQVTQQHNGW LMVDDAHGTG VIGEQGRGSC WLQKVKPELL VVTFGKGFV SGAAVLCSST VADYLLQFAR HLIYSTSMPP AQAQALRASL AVIRRDEGDA RREKLVSLIA RFRAGVQDLP FTLADSCSAI QPLIVGDNSR ALQLAEKLRQ QGCWVTAIRP PTVPAGTARL RLTLTAAHEM QDIDRLLEVL HGNG
Source	Mammalian cell
Target Names	bioF
Protein Names	Recommended name: 8-amino-7-oxononanoate synthase Short name= AONS EC= 2.3.1.47 Alternative name(s): 7-keto-8-amino-pelargonic acid synthase Short name= 7-KAP synthase Short name= KAPA synthase 8-amino-7- ketopelargo
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.