



# Recombinant Escherichia fergusonii Maltoporin (IamB)

<b>Product Code</b>	CSB-MP485019EOR
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7LL09
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	VDFHG YARSGIGWTG SGGEQQCFQT TGAQSKYRLG NECETYAELK LGQEVWKEGD KSFYFDTNVA YSVAQQNDWE ATSPAFREAN VQGKNLIDAL PGATIWAGKR FYQRHDVHMI DFYYWDISGP GAGLENIDLG FGKLSLAATR SSEAGGSSSF ASDNIYDYVN ETANDVFDVR LAQMEINPGG TLELGVDYGR ANLRDDYRLV DGASKDGWMF TAEHTQSMLK GFNKFVVQYA TDAMTSQGKG LNQGSGVAFD NEHFAYNINN NGHLLRILDH GAISLGDNWD MMYVGMQDI NWDNDNGTKW WTVGIRPMYK WTPIMSTLLE VGYDNVESQR TGDKNNQYKI TLAQQWQAGD SIWSRPAIRV FATYAKWDEK WGYDYNDNSK TNPNYGKAVP ANFEGGSFGR GDSDEWTFGA QMEIWW
<b>Source</b>	Mammalian cell
<b>Target Names</b>	IamB
<b>Protein Names</b>	Recommended name: Maltoporin Alternative name(s): Maltose-inducible porin
<b>Expression Region</b>	26-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.