



# Recombinant Escherichia fergusonii ATP-dependent Clp protease ATP-binding subunit ClpX (clpX)

<b>Product Code</b>	CSB-MP485030EOR
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7LME1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTDKRKDGSG KLLYCSFCGK SQHEVRKLIAGPSVYICDEC VDLNCNDIIRE EIKEVAPHRE RSALPTPHEI RHHLDDYVIG QEQAKKVLAV AVYNHYKRLR NGDTSNGVEL GKSNILLIGP TGSGKTLLE TLARLLDVPF TMADATTLTE AGYVGEDVEN IIQKLLQKCD YDVQKAQGI VYIDEIDKIS RKSDNPSITR DVS GEGVQQA LLKLI EGTVA AVPPQGG RKH PQQEFLQVDT SKILFICGGA FAGLDKVISH RVETGSGIGF GATVKAKSDK ASEGELLAQV EPEDLIKFG L IPEFIGRLPV VATLNELSEE ALIQILKEPK NALTKQYQAL FNLEGVDLEF RDEALDAIAK KAMARKTGAR GLRSIVEAAL LDTMYDLPSM EDVEKVVIDE SVIEGQSKPL LIYGKPEAQQ ASGE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	clpX
<b>Protein Names</b>	Recommended name: ATP-dependent Clp protease ATP-binding subunit ClpX
<b>Expression Region</b>	1-424
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.