



Recombinant Escherichia fergusonii Ribosomal protein S12 methylthiotransferase RimO (rimO)

| | |
|--------------------------|--|
| Product Code | CSB-YP485035EOR |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B7LMZ9 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSKVTHQPKI GFVSLGCPKN LVDSERILTE LRTEGYDVVP SYDDADMVIV NTCGFIDSAV QESLEAIGEA LNENGKVIVT GCLGAKEDQI REVHPKVLEI TGPHSYEQVL EHVHHYVPKP KHNPFSLVLP EQGVKLTPRH YAYLKISEGC NHRCTFCIIP SMRGDLVSRP IGDVLSEAKR LVDAGVKEIL VISQDTSAYG VDVKHRTGFY NGEPVK TSMV SLCEQLSKLG IWTRLHYVYP YPHVDDVIPL MAEGKILPYL DIPLQHASPR ILKLMKRP GS VDRQLARIKQ WREICPELTL RSTFIVGFPG ETEEDFQMLL DFLKEARLDR VGCFKYS PVE GADANALPDQ VPEEVKEERW NRFMQLQQQI SAERLQEKVG KEILVIIDEV DDEGAIGRSM ADAPEIDGAV YLNGESNVKP GDIVRVKVEH ADEYDLWGSR V |
| Source | Yeast |
| Target Names | rimO |
| Protein Names | Recommended name: Ribosomal protein S12 methylthiotransferase RimO Short name= S12 MTTase Short name= S12 methylthiotransferase EC= 2.-.-.- Alternative name(s): Ribosome maturation factor RimO |
| Expression Region | 1-441 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |