



# Recombinant Escherichia fergusonii Pyridoxine/pyridoxamine 5'-phosphate oxidase (pdxH)

<b>Product Code</b>	CSB-EP485052EOR-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7LQN2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSDNDELQQI AHLRREYTKG GLRRRDLPAD PLTLFERWLS QACEAKLADP TAMVVATVDE NGQPYQRIVL LKHYDEKGMV FYTNLGSRKA HQIESNPRVS LLFPWHTLER QVMVTGKAER LSTLEVVKYF HSRPRDSQIG AWWSKQSSRI SARGILESKF LELKQKFQQG EVPLPSFWGG FRVSMEQIEF WQGGEHRLHD RFLYQRENDK WKIDRLAP
<b>Source</b>	E.coli
<b>Target Names</b>	pdxH
<b>Protein Names</b>	Recommended name: Pyridoxine/pyridoxamine 5'-phosphate oxidase EC= 1.4.3.5 Alternative name(s): PNP/PMP oxidase Short name= PNPOx Pyridoxal 5'-phosphate synthase
<b>Expression Region</b>	1-218
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.