



Recombinant Escherichia fergusonii Protein CsiD (csiD)

Product Code	CSB-MP485069EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LT11
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MNALTAVQNN AVDSGQDYSYSG FTLIPSAQSP RLELTFTEQ TTKQFLEQVA EWPVQALEYK SFLRFRVGKI LDDLCANQLQ PLLLKTLLNR AEGALLINAV GVDDVAQADE MVKLATAVAH LIGRSNFDAM SGQYYARFVV KNVDNSDSYL RQPHRVMELH NDGTYYVEEIT DYVLMMKIDE QNMQGGNSLL LHLDDWEHLE HYFRHPLARR PMRFAAPPSK NVSKDVFHPV FDVDQQGRP V MRYIDQFVQP KDFEEGVWLS ELSDAIETSK GILSVPVPG KFLINLFW LHGRDRFTAH PDLRRELMRQ RGYFAYATHH YQTHQ
Source	Mammalian cell
Target Names	csiD
Protein Names	Recommended name: Protein CsiD
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.