



Recombinant Escherichia coli O8 Proline--tRNA ligase (proS)

Product Code	CSB-MP485139EEO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7M1Z9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O8 (strain IA11)
Purity	>85% (SDS-PAGE)
Sequence	MRTSQYLLSTLKETPADAEVISHQLMLRAGMIRKSLASGLYTWLPTGVRVLKKV ENIVREE MNNAGAIEVSMPVVQPADLWQESGRWEQYGPELLRFVDRGERPFVLGP THE EVITDLIRN ELSSYKQLPLNFYQIQTKFRDEVPRFRGVMRSREFLMKDAYSFHTSQESLQET YDAMYAA YSKIFSRMGLDFRAVQADTGSIGGSASHEFQVLAQSGEDDVVFSDTSDYAANI ELAEAIA PKEPRAAATQEMTLVDTPNAKTIAELVEQFNLPIEKTVKLLVKAVEGSSFPLVA LLVRG DHELNEVKAELPQVASPLTFATEEEIRAVVKAGPGSLGPVNMPPIPVIDRTVA AMSDFA AGANIDGKHYFGINWDRDVATPEVADIRNVVAGDPSPDGQGTLLIKRGIEVGH I FQLGTK YSEALKASVQGEDGRNQILTMGCYGIGVTRVVAAAIEQNYDERGIVWPDAIAP FQVAILP MNMHKSFRVQELAELKLYSELRAQGIEVLLDDRKERPGVMFADMELIGIPHTIVL GDRNLD NDDIEYKYRRNGEKQLIKTGDIVEYLVKQIKG
Source	Mammalian cell
Target Names	proS
Protein Names	Recommended name: Proline--tRNA ligase EC= 6.1.1.15 Alternative name(s): Prolyl-tRNA synthetase Short name= ProRS
Expression Region	1-572
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.