



# Recombinant Escherichia coli O45:K1 (Dimethylallyl)adenosine tRNA methylthiotransferase MiaB

<b>Product Code</b>	CSB-MP485243EOK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7MFS7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O45:K1 (strain S88 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTKKLHIKTW GCQMNEYDSS KMADLLDATH GYQLTDVAEE ADVLLLNTCS IREKAQEKVF HQLGRWKLK EKNPDLIGV GGCVASQEGE HIRQRAHYVD IIFGPQTLHR LPEMINSVRG DRSPVVDISF PEIEKFDRLP EPRAEGPTAF VSIMEGCKNY CTYCVVPYTR GEEVSRPSDD ILFEIAQLAA QGVREVNLLG QNVNAWRGEN YDGTGTSFAD LLRLVAAIDG IDRIRFTTSH PIEFTDDIIE VYRDTPELVS FLHLPVQSGS DRILNLMGRT HTALEYKAI RKLRAARPDI QISSDFIVGF PGETTDDFEK TMKLIADVNF DMSYSFIFSA RPGTPAADMV DDVPEEEKKQ RLYILQERIN QQAMAWSRRM LGTTQRILVE GTSRKSIMEL SGRTENNRVV NFEGTPDMIG KFVDVEITDV YPNSLRGKVV RTEDEMGLRV AETPESVIAR TRKENDLGVG YYQP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	miaB
<b>Protein Names</b>	Recommended name: (Dimethylallyl)adenosine tRNA methylthiotransferase MiaB EC= 2.-.-.- Alternative name(s): tRNA-i(6)A37 methylthiotransferase
<b>Expression Region</b>	1-474
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.