



Recombinant Escherichia coli O45:K1 Succinyl-diaminopimelate desuccinylase (dapE)

Product Code	CSB-BP485260EOK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7MHW4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O45:K1 (strain S88 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSCPVIELTQ QLIRRPSLSP DDAGCQALLI ERLQAIGFTV ERMDFADTQN FWAWRGQGET LAFAGHTDVV PPGDADRWIN PPFPTIRDG MLFGRGAADM KGSLAAMVVA AERFVAQHPN HTGRLAFLIT SDEEASAHNG TVKVVEALMA RNERLDYCLV GEPSSIEVVG DVVKNRRGS LTCNLTIHGV QGHVAYPHLA DNPVHRAAPF INELVAIEWD QGNEFFPATS MQIANIQAGT GSNNVIPGEL FVQFNFRFST ELTDEMIKAQ VLALLEKHQL RYTVDDWWLSG QPFLTARGKL VDAVVNAVEH YNEIKPQLLT TGGTSDGRFI ARMGAQVVEL GPVNATIHKI NECVNAADLQ LLARMYQRIM EQLVA
Source	Baculovirus
Target Names	dapE
Protein Names	Recommended name: Succinyl-diaminopimelate desuccinylase Short name= SDAP desuccinylase EC= 3.5.1.18 Alternative name(s): N-succinyl-L-lysine desuccinylase
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.