



# Recombinant Escherichia coli O81 N-methyl-L-tryptophan oxidase (solA)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP485323EOP   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B7MTJ0  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O81 (strain ED1a)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MKYDLIIIGS GSVGAAAGYY ATRAGLNVLM TDAHMPPHQH GSHHG DTRLI<br>RHAYGEGEKY VPLVLRQML WDELSRHNE DPFVRSQVI NLGPADSAFL<br>ANVAHSAEQW QLNVEKLDAQ GIMARWPEIR VPDNYIGLFE TDSGFLRSEL<br>AIKTWQLAK EAGCAQLFNC PVTAIRHDDD GVTIETADGV YQAKKAIVCA<br>GTWVKDLLPE LPVQPVRKVF AWYQADGRYS VKNKFFAFTG ELPNGDQYYG<br>FPAENDALKI GKHNGGQVIH SADERVPFAE VVSDGSEAFP FLRNVLPGIG<br>CCLYGAACTY DNSPDEDFII DTLPGHDNTL LITGLSGHGF KFASVLGEIA<br>ADFAQDKKSD FDLTPFRLSR FQ |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | solA  |
| <b>Protein Names</b>     | Recommended name: N-methyl-L-tryptophan oxidase Short name= MTOX EC= 1.5.3.-  |
| <b>Expression Region</b> | 1-372   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                 |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |