



Recombinant Escherichia coli O17:K52:H18 Peptide chain release factor 2 (prfB)

Product Code	CSB-MP485413EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7N7D5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MFEINPVNNR IQDLTERS DV LRGYLDYDAK KERLEE VNAE LEQPDVWNEP ERAQALGKER SSLEAVVDTL DQMKQGLEDV SGLLELAVEA DDEETFNEAV AELDALEEKL AQLEFRRMFS GEYDSADCYL DIQAGSGGTE AQDWASMLER MYLRWAESRG FKTEIIEESE GEVAGIKSVT IKISGDYAYG WLRTETGVHR LVRKSPFDSG GRRHTSFSSA FVYPEVDDDI DIEINPADLR IDVYRASGAG GQHVNRTESA VRITHIPTGI VTQCQNDRSQ HKNKDQAMKQ MKAKLYELEM QKKNAEKQAM EDNKSDIGWG SQIRSYVLDD SRIKDLRTGV ETRNTQAVLD GSLDQFIEAS LKAGL
Source	Mammalian cell
Target Names	prfB
Protein Names	Recommended name: Peptide chain release factor 2 Short name= RF-2
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.