



Recombinant Escherichia coli O17:K52:H18 Deoxyguanosinetriphosphate triphosphohydrolase (dgt)

Product Code	CSB-BP485423EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7N829
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MAQIDFRKKI NWHRRYRSPQ GVKTEHEILR IFESDRGR II NSPAIRRLQQ KTQVFPLERN AAVRTRLTHS MEVQQVGRIYI AKEILSRLKE LKLEAYGLD ELTGPFESIV EMSCLMHDIG NPPFGHFGEA AINDWFRQRL HPEDAESQPL TDDRCSVAAL RLRDGEPLN ELRRKIRQDL CHFEGNAQGI RLVHTLMRMN LTWAQVGGIL KYTRPAWWRG ETPETHHYLM KKPGYYLSEE AYIARLRKEL NLALYSRFPL TWIMEAADDI SYCVADLEDA VEKRIFTVEQ LYHHLHEAWG QHEKGSLSFSL VVENAWEKSR SNSLSRSTED QFFMYLRVNT LNKLVYAAQ RFIDNLPAIF AGTFNHALLE DASECSDLLK LYKNVAVKHV FSHPDVEQLE LQGYRVISGL LEIYRPLLSL SLSDFTELVE KERVKRFPPIE SRLFNKLSTR HRLAYVEAVS KLPSDSPEFP LWEYYYRCRL LQDYISGMTD LYAWDEYRRL MAVEQ
Source	Baculovirus
Target Names	dgt
Protein Names	Recommended name: Deoxyguanosinetriphosphate triphosphohydrolase Short name= dGTP triphosphohydrolase Short name= dGTPase EC= 3.1.5.1
Expression Region	1-505
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.