



# Recombinant Escherichia coli O17:K52:H18 1-deoxy-D-xylulose-5-phosphate synthase (dxs)

<b>Product Code</b>	CSB-MP485430EOG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7N8X3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSFDIAKYPT LALVDSTQEL RLLPKESLPK LCDELRRYLL DSVSRSSGHF ASGLGTVELT VALHYVYNTF FDQLIWDVGH QAYPHKILTG RRDKIGTIRQ KGGLHPFPWR GESEYDVLSV GHSSTSISAG IGIAVAAEKE GKNRRTVCVI GDGAITAGMA FEAMNHAGDI RPDMLVVLND NEMSIENVG ALNNHLAQLL SGKLYSSLRE GGKKVFSGVP PIKELLK RTE EHIKGMVVP G TLFEELGFNY IGPVDGHDVL GLITTLKNMR DLKGPQFLHI MTKKGRGYEP AEKDPITFHA VPKFDPSGC LPKSSGGLPS YSKIFGDWLC ETAAKDNKLM AITPAMREGS GMVEFSRKFP DRYFDVAIAE QHAVTFAAGL AIGGYKPVVA IYSTFLQRAY DQVLHDVAIQ KLPVLF AIDR AGIVGADGQT HQGAFDLSYL RCIPEMVIMT PSDENECRQM LYTG YHYNDG PSAVRYPRGN AVGVELTPLE KLPIGKGIVK RRGEKLAILN FG TLMPEAAK VAESLNATLV DMR FVKPLDE ALILEMAASH EALVTVEENA IMGGAGSGVN EVLMAHRKPV PVLNIGLPDF FIPQGTQEEM RAELSLDAAG MEAKIKAWLA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	dxs
<b>Protein Names</b>	Recommended name: 1-deoxy-D-xylulose-5-phosphate synthase EC= 2.2.1.7 Alternative name(s): 1-deoxyxylulose-5-phosphate synthase Short name= DXPS synthase Short name= DXPS
<b>Expression Region</b>	1-620
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.