



# Recombinant Escherichia coli O17:K52:H18 Octanoyltransferase (lipB)

<b>Product Code</b>	CSB-EP485434EOG-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7N9N6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MYQDKILVRQ LGLQPYEPIS QAMHEFTDTR DESTLDEIWL VEHYPVFTQG QAGKAEHILM PGDIPVIQSD RGGQVTYHGP GQQVMYVLLN LKRRKLGVRE LVTLLLEQTVV NTLAELGIEA HPRADAPGVY VGEKKICSLG LRIRRGCSFH GLALNVNMDL SPFLRINPCG YAGMEMAKIS QWKPEATTNN IAPRLLENIL ALLNNPDFEY ITA
<b>Source</b>	E.coli
<b>Target Names</b>	lipB
<b>Protein Names</b>	Recommended name: Octanoyltransferase EC= 2.3.1.181 Alternative name(s): Lipoate-protein ligase B Lipoyl/octanoyl transferase Octanoyl-[acyl-carrier- protein]-protein N-octanoyltransferase
<b>Expression Region</b>	1-213
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.