



Recombinant Escherichia coli O17:K52:H18 Phosphopentomutase (deoB)

Product Code	CSB-EP485486EOG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NH51
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKRAFIMVLD SFGIGATEDA ERFGDVGADT LGHIAEACAK GEADNGRKGP LNLPNLTRLG LAKAHEGSTG FIPAGMDGNA EVIGAYAWAH EMSSGKDTPS GHWEIAGVPV LFEWGYFSDH ENSFPQELLD KLVERANLPG YLGNCHSSGT VILDQLGEEH MKTGKPIFYT SADSVFQIAC HEETFGLDKL YELCEIAREE LTNGGYNIGR VIARPFIGDK AGNFQRTGNR HDLAVEPPAP TVLQKLVDK HGQVSVGKI ADIYANCGIT KVKATGLDA LFDATIKEMK EAGDNTIVFT NFVDFDSSWG HRRDVAGYAA GLELFDRRLP ELMSLLRDDD ILILTADHGC DPTWTGTDHT REHIPVLVYG PKVKPGSLGH RETFADIGQT LAKYFGTSDM EYGKAMF
Source	E.coli
Target Names	deoB
Protein Names	Recommended name: Phosphopentomutase EC= 5.4.2.7 Alternative name(s): Phosphodeoxyribomutase
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.