



Recombinant Escherichia coli O7:K1 Hydroxyethylthiazole kinase (thiM)

Product Code	CSB-BP485540EON
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NPM9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MQVDLLSSAQ SAHALHLFHK HSPLVHCMTN DVVQTFTANT LLALGASPAM VIETEEASQF AAIASALLIN VGTLTQPRAQ AMSAAVEQAT RSQTPWTLDP VAVGALDYRR RFCVELLSHK PTAIRGNASE IMALAGVANG GRGVDTTDAA ANAIPAAQTL ARETGAIVVV TGEVDYVTDG HRIVGIHGGD PLMTKVVG TG CALSAVVAAC CALPGDTLEN IASACHWMKQ AGERAVARSE GPGSFVPHFL DALWQLTQEV QA
Source	Baculovirus
Target Names	thiM
Protein Names	Recommended name: Hydroxyethylthiazole kinase EC= 2.7.1.50 Alternative name(s): 4-methyl-5-beta-hydroxyethylthiazole kinase Short name= TH kinase Short name= Thz kinase
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.