



# Recombinant Escherichia coli O7:K1 Chromosomal replication initiator protein DnaA (dnaA)

<b>Product Code</b>	CSB-MP485544EON
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7NR04
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSLSLWQQCL ARLQDELPAT EFSMWIRPLQ AELSDNTLAL YAPNRFVLDW VRDKYLNNIN GLLTSFCGAD APQLRFEVGT KPVTQTPQAA VTSNVAAPAAQ VAQTQPQRAA PSTRSGWDNV PAPAPTYRS NVNVKHTFDN FVEGKSNQLA RAAARQVADN PGGAYNPLFL YGGTGLGKTH LLHAVGNGIM ARKPNKVVY MHSERFVQDM VKALQNNIAIE EFKRYYSVD ALLIDDIQFF ANKERSQEEF FHTFNALLEG NQQIILTSR YPKEINGVED RLKSRFGWGL TVAIEPPELE TRVAILMKKA DENDIRLPGE VAFFIAKRLR SNVRELEGAL NRVIANANFT GRAITIDFVR EALRDLLALQ EKLVTIDNIQ KTVAEYKIK VADLLSKRRS RSVARPRQMA MALAKELTNH SLPEIGDAFG GRDHTTVLHA CRKIEQLREE SHDIKEDFSN LIRTLSS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	dnaA
<b>Protein Names</b>	Recommended name: Chromosomal replication initiator protein DnaA
<b>Expression Region</b>	1-467
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.