



Recombinant Escherichia coli O7:K1 Nitric oxide reductase FIRd-NAD (+) reductase

Product Code	CSB-MP485556EON
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NSJ1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSNSIVIIGS GFAARQLVKN IRKQDASIPL TLIAADSMDE YNKPDLSHVI SQGQRADDLT RQTAGEFAEQ FNLRLFPHW VTDIDAEAHV VKSQNNQWQY DKLVLATGAS AFVPPVPGRE LMLTLNSQQE YRACETQLRD ARRVLIVGGG LIGIELAMDF CRAGKAVTLI DNAASILASL MLPEVSSRLQ HRLTEMGVHL LLKSQLQGLE KTDGILATL DRQRSIEVDA VIAATGLRPE TALARRAGLT INRGVCVDSY LQTSNADIYA LGDCAEINGQ VLPFLQPIQL SAMVLAKNLL GNNTPLKLPALMLVKIKTPEL PLHLAGETQR QDLRWQINTE RQGMVARGVD DGDQLRAFVV SEDRMKEAFG LLKTLPM
Source	Mammalian cell
Target Names	norW
Protein Names	Recommended name: Nitric oxide reductase FIRd-NAD(+) reductase EC= 1.18.1.- Alternative name(s): Flavorubredoxin reductase Short name= FIRd-reductase Short name= FlavoRb reductase
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.