



# Recombinant Escherichia coli O7:K1 Anhydro-N-acetylmuramic acid kinase (anmK)

<b>Product Code</b>	CSB-YP485567EON
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7NTZ2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKSGRFIVGM SGTSLDGVVDV VLATIDEHRV AQLASLSWPI PVSLKQAVLD ICQQQQLTLS QFGQLDTQLG RLFADAVNAL LKEQNLQARD IVAIGCHGQT VWHEPTGVAP HTLQIGDNNQ IAARTGITVV GDFRRRDIAL GGQGAPLVPA FHHALLAHPT ERRMVLNIGG IANLSLLIPG QPVG GYDTGP GNMLMDAWIW RQAGKPYDKD AEWARAGKVI LPLLQNMLSD PYFSQPAPKS TGREYFN YGW LERHLRHFPV VDP RDVQATL AELTAVTISE QVLLSGGCER LMVCGGGSRN PLLMARLAAL LPGTEVTTTD AVGISGDDME ALAFAWLAWR TLAGLPGNLP SVTGASQETV LGAIFPANP
<b>Source</b>	Yeast
<b>Target Names</b>	anmK
<b>Protein Names</b>	Recommended name: Anhydro-N-acetylmuramic acid kinase EC= 2.7.1.170 Alternative name(s): AnhMurNAc kinase
<b>Expression Region</b>	1-369
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.