



Recombinant Escherichia coli O127:H6 4-hydroxy-3-methylbut-2-en-1-yl diphosphate synthase (ispG)

Product Code	CSB-BP485594EOB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7UGW1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Purity	>85% (SDS-PAGE)
Sequence	MHNQAPIQRR KSTRIYVGNV PIGDGAPIAV QSMTNTRTTD VEATVNQIKA LERVGADIVR VSVPTMDAAE AFKLIKQQVN VPLVADIHFD YRIALKVAEY GVDCLRINPG NIGNEERIRM VVDCARDKNI PIRIGVNAGS LEKDLQEKYG EPTPQALLES AMRHVDHLDR LNFDQFKVSV KASDVFLAVE SYRLAKQID QPLHLGITEA GGARSGAVKS AIGLGLLLSE GIGDTRLVSL AADPVVEEIKV GFDILKSLRI RSRGINFIAC PTCRQEFDV IGTVNALEQR LEDIITPMDV SIIGCVVNGP GEALVSTLGV TGGNKKSGLY EDGVRKDRLD NNDMIDQLEA RIRAKASQLD EARRIDVQQV EK
Source	Baculovirus
Target Names	ispG
Protein Names	Recommended name: 4-hydroxy-3-methylbut-2-en-1-yl diphosphate synthase EC= 1.17.7.1 Alternative name(s): 1-hydroxy-2-methyl-2-(E)-butenyl 4-diphosphate synthase
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.