



# Recombinant Escherichia coli O127:H6 Glycerol kinase (glpK)

<b>Product Code</b>	CSB-MP485649EOB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7UNP6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTEKKYIVAL DQGTSSRAV VMDHDANIIS VSQREFEQIY PKPGWVEHDP MEIWATQSST LVEVLAKADI SSDQIAAIGI TNQRETTIVW EKETGKPIYN AIVWQCRRTA EICEHLKRDG MEEYIRNNTG LVIDPYFSGT KVKWILDHVE GSRERARRGE LLFGTVDTWL IWKMTQGRVH VTDYTNASRT MLFNIHTLDW DDKMLEVLDI PREMLPEVRR SSEVYGQTNI GGKGGTRIPi SGIAGDQQAA LFGQLCVKEG MAKNTYGTGC FMLMNTGEKA VKSENGLLTT IACGPTGEVN YALEGAVFMA GASIQWLRDE MKLINDAYDS EYFATKVQNT NGVYVVAFT GLGAPYWDPY ARGAI FGLTR GVNANHIIRA TLESIAYQTR DVLEAMQADS GIRLHALRVD GGAVANNFLM QFQSDILGTR VERPEVREVT ALGAAYLAGL AVGFWQNLDE LQEKAVIERE FRPGIETTER NYRYAGWKKA VKRAMAWEEH DE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	glpK
<b>Protein Names</b>	Recommended name: Glycerol kinase EC= 2.7.1.30 Alternative name(s): ATP:glycerol 3-phosphotransferase Glycerokinase Short name= GK
<b>Expression Region</b>	1-502
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.