



Recombinant Escherichia fergusonii Erythronate-4-phosphate dehydrogenase (pdxB)

Product Code	CSB-EP485858EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LLF0
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MKILVDENMP YARELFSRLG EVKAVPGRPI PVAQLADADA LMVRSVTKVN ESLLAGKPIK FVGTATAGTD HVDEAWLKQA GIGFSAAPGC NAIAVVEYVF SLLMLAERD GFSLHERTVG IVGVGNVGRR LQARLEALGI KTLLCDPPRA DRGDEGDFRS LDELVQHADI LTFHTPLFKD GPYKTLHLAD EKLIRSLKPG AILINACRGA VVDNTALLTC LNEGQKLSVV LDVWEGEPEL NVELLKKVDI GTPHIAGYTL EGKARGTTQV FEAYSKFIGH EQHVALDTLL PAPEFGRITL HGLLDQPTLK RLVHLVYDVR RDDAPLRKVA GIPGEFDKLR KNYLERREWS SLYVICDDAS AASLLCKLGF NAVHHPAR
Source	E.coli
Target Names	pdxB
Protein Names	Recommended name: Erythronate-4-phosphate dehydrogenase EC= 1.1.1.290
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.