



Recombinant Escherichia fergusonii Membrane-bound lytic murein transglycosylase C (mltC)

Product Code	CSB-BP485881EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LPT3
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	≥85% (SDS-PAGE)
Sequence	CSTT KKGDTYNEAW VKDTNGFDIL MGQFAHNIEN IWGFKEVVIA GPKDYVKYTD QYQTRSHINF DDGTITIETI AGTEPAAHLR RAIKTLMLG DDPSSVDLYS DVDDITISKE PFLYGQVVDN TGQPIRWEGR ASNFADYLLK NRLKSRSNGL RIIYSVTINM VPNHLDKRAH KYLGMVRQSS RKYGVDESLE LAIMQTESSF NPYAVSRSDA LGLMQVVQHT AGKDVFRSQG KSGTPSRNFL FDPASNIDTG TAYLAMLNNV YLGGIDNPTS RRYAVITAYN GGAGSVLRVF SNDKIQAANI INTMTPGDVY QTLTTRHPSA ESRRYLYKVN TAQKSYRRR
Source	Baculovirus
Target Names	mltC
Protein Names	Recommended name: Membrane-bound lytic murein transglycosylase C EC=3.2.1.- Alternative name(s): Murein hydrolase C
Expression Region	17-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.