



Recombinant Escherichia fergusonii N-succinylglutamate 5-semialdehyde dehydrogenase (astD)

Product Code	CSB-EP485882EOR-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LQ46
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MTLWINGDWV TGQGALRVKC NPVSGELLWQ GNDADAAQVG QACRAARA AF PRWARLSFGD RQVRVERFAG LLESNKVELT AIIARETGKP RWEAATEVTA MINKIAISIK AYHVRTGEQR SEMPDGAASL RHRPHGV LAV FGPYNFPGHL PNGHVPALL AGNTVIFKPS ELTPWSGEAV MRLWQQAGLP SGVLNLVQGG RETGQALSAL EDLDG LLFTG SANTGYQLHR QLSGQPEKIL ALEMGGNNPL IIDDVADIDA AVHLTIQSAF VTAGQRCTCA RRLLLKSGAQ GDAFLARLVA VSQRLTPGNW DDEPQPFIGG LISEQAAHQV VTAWQELEAM GGRTQLAPRL LRAGTSL LTP GIVEMTG VAG VPDEEVFGPL LRVWRYDTFD EAIMANNTR FGLSCGLVSS RRDKFEQLLL EARAGIVNWN KPLTGAASTA PFGVGASGN HRPSAWYAAD YCAWPMASLE SDSLTLPTTL NPGLDFSEEV ER
Source	E.coli
Target Names	astD
Protein Names	Recommended name: N-succinylglutamate 5-semialdehyde dehydrogenase EC= 1.2.1.71 Alternative name(s): Succinylglutamic semialdehyde dehydrogenase Short name= SGSD
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.