



Recombinant Escherichia fergusonii N-methyl-L-tryptophan oxidase (solA)

Product Code	CSB-YP485907EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LT77
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	≥85% (SDS-PAGE)
Sequence	MKYDLIIIGS GSVGAAAGYY ATRAGLNVLM IDAHMPPHQQ GSHHGDTRLI RHAYGEGEKY VPLVLRQQQL WDELSTHNEE PVFVRSGVIN LGPENSEFLA NVAHSAEQWQ LPVERLDATA ITSRWPEIRV PENYIGLFES ESGFLRSELA IETWIRLAKE AGCAQLFNCP VSAIRHEGDG VTVVTEDEGEY SARKVLISAG TWVKSLLEL PVQPVKVFV WYQADGRYSI KKNKFPFTGE LPNGDQYYGF PAENDALKIG KHNGGQSISN AAERLPFAEV ASDGSEAFPFP LRSVLPGIGC CLYGASCTYD NSPDEDFIID TLPGHDNTLV ITGLSGHGFK FASVLGEIAA DFAQGKPAAF DLTPFKLSRF K
Source	Yeast
Target Names	solA
Protein Names	Recommended name: N-methyl-L-tryptophan oxidase Short name= MTOX EC= 1.5.3.-
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.