



Recombinant Escherichia fergusonii Bis (5'-nucleosyl)-tetrphosphatase, symmetrical

| | |
|--------------------------|---|
| Product Code | CSB-MP485925EOR |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B7LVU3 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MATYLIGDVH GCYDELIALL HKVEFTPGKD TLWLTGDLVA RGPGLDVLRL YVKSLGDSVR LVLGNHDLHL LAVFAGISRN KPKDRLTPLL EAPDADELLN WLRRQPLLQI DEEKKLVMAH AGITPQWDLQ TAKECARDVE AVLSSDSYPF FLDAMYGDMP NNWSPELRGL GRLRFITNAF TRMRFCFPNG QLDMYSKESP EEAPAPLKPW FAIPGPVAEE YSIAFGHWAS LEGKGTPEGI YALDTGCCWG GTLTCLRWED KQYFVQPSNR HKDLGEGEAV AS |
| Source | Mammalian cell |
| Target Names | apaH |
| Protein Names | Recommended name: Bis(5'-nucleosyl)-tetrphosphatase, symmetrical EC= 3.6.1.41 Alternative name(s): Ap4A hydrolase Diadenosine 5',5'''-P1,P4- tetrphosphate pyrophosphohydrolase Diadenosine tetrphosphatase |
| Expression Region | 1-282 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |