



Recombinant Escherichia fergusonii Glutamate-1-semialdehyde 2,1-aminomutase (hemL)

Product Code	CSB-BP485931EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LWB5
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MSKSENLISA ARELIPGGVN SPVRAFTGVG GTPLFIEKAD GAYLYDVDGK AYIDYVGSWG PMVLGHNHPA IRNAVIEAAE RGLSFGAPTE MEVKMAQLVT ELVPTMDMVR MVNSGTEATM SAIRLARGFT GRDKIIFEG CYHGHADCLL VKAGSGALTL GQPNSPGVPA DFAKHTLTCT YNDLASVRAA FEQYPQEIAC IIVEPVAGNM NCVPLPEFL PGLRALCDEF GALLIIDEVM TGFRVALAGA QDYDVPDL TCLGKIIGGG MPVGAFGGRR DVMDALAPIG PVYQAGTLSG NPIAMAAGFA CLNEVAQPGV HETLDELTTTR LAEGLLEAAE EAGIPLVVNH VGGMFGIFFT DAESVTCYQD VMACDVERFK RFFHMMLDEG VYLAPSAFEA GFMSVAHSME DINNTIDAAR RVFAKL
Source	Baculovirus
Target Names	hemL
Protein Names	Recommended name: Glutamate-1-semialdehyde 2,1-aminomutase Short name= GSA EC= 5.4.3.8 Alternative name(s): Glutamate-1-semialdehyde aminotransferase Short name= GSA-AT
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.