



Recombinant Mouse Transcription factor jun-D (Jund)

Product Code	CSB-EP011974MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15066
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	METPFYGEEA LSGLAAGASS VAGATGAPGG GGFAPPGRAF PGAPPTSSML KKDALTLSLA EQGAAGLKPG SATAPSALRP DGAPDGLLAS PDLGLLKLAS PELERLIQS NGLVTTTPTS TQFLYPKVAA SEEQEFAEGF VKALEDLHKQ SQLGAATAAT SGAPAPPAPA DLAATPGATE TPVYANLSSF AGGAGPPGGA ATVAFAAEPV PFPPPGALG PPPPPHPPRL AALKDEPQTV PDVPSFGDSP PLSPIDMDTQ ERIKAERKRL RNRIAASKCR KRKLERISRL EEKVKTLKSQ NTELASTASL LREQVAQLKQ KVLSHVNSGC QLLPQHQPVA Y
Source	E.coli
Target Names	Jund
Protein Names	Recommended name: Transcription factor jun-D
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The protein encoded by this intronless gene is a member of the JUN family, and a functional component of the AP1 transcription factor complex. It has been proposed to protect cells from p53-dependent senescence and apoptosis. Alternate translation initiation site usage results in the production of different isoforms.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.