



Recombinant Escherichia coli O8 4-hydroxythreonine-4-phosphate dehydrogenase (pdxA)

Product Code	CSB-MP485957E00
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7M0E9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O8 (strain IA11)
Purity	>85% (SDS-PAGE)
Sequence	MVKTQRVVIT PGEPAGIGPD LVVQLAQREW PVELVVCADA TLLTDRAAML GLPLTLRPYS PNSPAQPQT GTLTLLPVAL RESVTAGQLA IENGHYVVET LARACDGCLN GEFAALITGP VHKGVINDAG IPFTGHTEFF EERSQAKKVV MMLATEELRV ALATTHLPLR DIADAITPAL LHEVAILHH DLRTKFGIAE PRILVCGLNP HAGEGGHMG T EEIDTIIPVL DELRAQGMKL NGPLPADTLF QPKYLDNADA VLAMYHDQGL PVLKYQGFR GVNITLGLPF IRTSVDHGTA LELAGRGEAD VGSFITALNL AIKMIVNTQ
Source	Mammalian cell
Target Names	pdxA
Protein Names	Recommended name: 4-hydroxythreonine-4-phosphate dehydrogenase EC=1.1.1.262 Alternative name(s): 4-(phosphohydroxy)-L-threonine dehydrogenase
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.