



Recombinant Escherichia coli O45:K1 Fe/S biogenesis protein NfuA (nfuA)

Product Code	CSB-MP486066EOK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7MDP0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O45:K1 (strain S88 / ExPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MIRISDAAQA HFAKLLANQE EGTQIRVFVI NPGTPNAECG VSYCPPDAVE ATDTALKFDL LTAYVDELSA PYLEDAEIDF VTDQLGSQLT LKAPNAKMRK VADDAPLMER VEYMLQSQIN PQLAGHGGRV SLMEITEDGY AILQFGGGCN GCSMVDVTLK EGIEKQLLNE FPELKGVRDL TEHQRGEHSY Y
Source	Mammalian cell
Target Names	nfuA
Protein Names	Recommended name: Fe/S biogenesis protein NfuA
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.