



Recombinant Escherichia coli O45:K1 D-tagatose-1,6-bisphosphate aldolase subunit GatY (gatY)

| | |
|--------------------------|---|
| Product Code | CSB-EP486071EOK-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B7MEF1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O45:K1 (strain S88 / ExPEC) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MYVVSTKQML NNARRGGYAV PAFNIHNLET MQVVTETAAS MHAPVIIAGT PGTFTHAGTE NLMALVSAMA KQYHHPLAIH LDHHTKFDDI AQKVRSGVRS VMIDASHLPF AQNISRVKEV VDFCHRFDVS VEAELGQLGG QEDDVQVNEA DAFYTNPVQA REFAEATGID SLAVAIGTAH GMYASAPALD FSRLNIRQW VNLPLVLHGA SGLSTKDIQQ TIKLGICKIN VATELKNAPS QALKNYLTEY PEATDPRDYL QSAKSAMRDV VSKVIADCGC EGRA |
| Source | E.coli |
| Target Names | gatY |
| Protein Names | Recommended name: D-tagatose-1,6-bisphosphate aldolase subunit GatY Short name= TBPA Short name= TagBP aldolase EC= 4.1.2.40 Alternative name(s): D-tagatose-bisphosphate aldolase class II Tagatose-bisphosphate aldolase |
| Expression Region | 1-284 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |