



Recombinant Escherichia coli O17:K52:H18 HTH-type transcriptional regulator BetI (betI)

Product Code	CSB-MP486271EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7N8L5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MPKLGMQSIR RRQLIDATLE AINEVGMHDA TIAQIARRAG VSTGIISHYF RDKNGLLEAT MRDITSQLRD AVLNRLHALP QGSAEQRLQA IVDGNFDETQ VSNAAMKAWL AFWASSMHQP MLYRLQQVSS RRLLSNLVSE FHRELPRQQA QEAGYGLAAL IDGLWLRAAL SGKPLDKPLA HSLTRHFITQ HLPTD
Source	Mammalian cell
Target Names	betI
Protein Names	Recommended name: HTH-type transcriptional regulator BetI
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.