



Recombinant Escherichia coli O17:K52:H18 Histidinol-phosphate aminotransferase (hisC)

Product Code	CSB-MP486292EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NC61
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSTVTITDLA RENVRLTPY QSARRLGGNG DVWLNANEYP TAVEFQLTQQ TLNRYPECQP KAVIENYAQY AGVKPEQVLV SRGADEGIEL LIRAFCEPGK DAILYYPPTY GMYSVSAETI GVECRTVPTL DNWQLDLQGI SDKLDGVKVV YVCSNNPTG QLINPQDFRI LLELTRGKAI VVADEAYIEF CPQASLAGWL TEYPHLAILR TLSKAFALAG LRCGFTLANE EVINLLMKVI APYPLSTPVA DIAAQALSPQ GIVAMRERVA QIIAEREYLI AALKEISCVE QVFDSETNYI LARFKASSAV FKSLWDQGII LRDQNKQPSL SGCLRITVGT REESQRVIDA LRAEQV
Source	Mammalian cell
Target Names	hisC
Protein Names	Recommended name: Histidinol-phosphate aminotransferase EC= 2.6.1.9 Alternative name(s): Imidazole acetol-phosphate transaminase
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.