



Recombinant Escherichia coli O17:K52:H18 Putative N-acetylmannosamine-6-phosphate 2- epimerase (nanE)

Product Code	CSB-YP486299EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NDK1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MSLLAQLDQK IAANGGLIVS CQPVPDSPLD KPEIVAAMAL AAEQAGAVAI RIEGVANLQA TRAVVSVPII GIVKRDLEDS PVRITAYIED VDALAQAGAD IIAIDGTDRP RPVPVETLLA RIHHHGLLAM TDCSTPEDGL ACQKLGAEII GTTLSGYTTP ETPEEPDLAL VKTLSDAGCR VIAEGRYNTP AQAADAMRHG AWAVTVGSAI TRLEHICQWY NTAMKKAVL
Source	Yeast
Target Names	nanE
Protein Names	Recommended name: Putative N-acetylmannosamine-6-phosphate 2-epimerase EC= 5.1.3.9 Alternative name(s): ManNAc-6-P epimerase
Expression Region	1-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.