



# Recombinant Escherichia coli O17:K52:H18 Glucose-1-phosphate adenylyltransferase (glgC)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-YP486307EOG  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | B7NE40   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MVSLEKNDHL MLARQLPLKS VALILAGGRG TRLKDLTNKR AKPAVHFGGK<br>FRIIDFALSN CINSGIRRMG VITQYQSHTL VQHIQRGWSF FNEEMNEFVD<br>LLPAQQRMKG ENWYRGTADA VTQNLDIIRR YKAEYVVILA GDHIYKQDYS<br>RMLIDHVEKG ARCTVACMPV PIEEASAFGV MAVDENDKII EFVEKPANPP<br>SMPNDPSKSL ASMGYVFDA DYLYELLEED DRDENSSSHDF GKDLIPKITE<br>AGLAYAHPFP LSCVQSDPDA EPYWRDVGTL EAYWKANLDL ASVVPEDMY<br>DRNWPIRTYN ESLPPAKFVQ DRSGSHGRTL NSLVSGGCVI SGSVVVQSVL<br>FSRVRVNSFC NIDSAVLLPE VVVGRCRLR RCVIDRACVI PEGMVIGENA<br>EEDARRFYRS EEGIVLVTRE MLRKLGHKQE R |
| <b>Source</b>            | Yeast  |
| <b>Target Names</b>      | glgC   |
| <b>Protein Names</b>     | Recommended name: Glucose-1-phosphate adenylyltransferase EC= 2.7.7.27<br>Alternative name(s): ADP-glucose pyrophosphorylase Short name= ADPGlc<br>PPase ADP-glucose synthase  |
| <b>Expression Region</b> | 1-431  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |