



Recombinant Escherichia coli O7:K1 2-C-methyl-D-erythritol 4-phosphate cytidyltransferase (ispD)

Product Code	CSB-MP486407EON
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NT91
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MATTHLDVCA VVPAAGFGRR MQTECPKQYL SIGNQTILEH SVHALLAHPR VKHVVIAISP GDSRFAQQLPL ANHPQITVVD GGEERADSVL AGLKAAGDAQ WVLVHDAARP CLHQDDLARL LALSETSRTG GILAAPVRDT MKRAEPGKNA IAHTVDRNGL WHALTPQFFP RELLHDCLTR ALNEGATITD EASALEYCGF HPQLVEGRAD NIKVTRPEDL ALAEFYLTRT IHQENT
Source	Mammalian cell
Target Names	ispD
Protein Names	Recommended name: 2-C-methyl-D-erythritol 4-phosphate cytidyltransferase EC= 2.7.7.60 Alternative name(s): 4-diphosphocytidyl-2C-methyl-D-erythritol synthase MEP cytidyltransferase Short name= MCT
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.