



Recombinant Escherichia coli O7:K1 UPF0061 protein ydiU (ydiU)

Product Code	CSB-BP486412EON
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NTS5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MTLSFITRWC DELPETYTAL SPTPLNKARL IWHNAELANT LSIPSSLFKN GAGVWGGETL LPGMSPLAQV YSGHQFGVWA GQLGDGRGIL LGEQQLADGT TMDWHLKGAG LTPYSRMGDG RAVLRSTIRE SLASEAMHYL GIPTRALSI VTSDSPVYRE TVEPGAMLMR VAPSHLRFGH FEHFYRREP EKVRQLADFA IRHYWSHLAD DEDKYRLWFS DVVARTASLI AQWQTVGFAH GVMNTDNMSL LGLTLDYGPF GFLDDYEPGF ICNHSDHQGR YSFDNQPAVA LWNLQRLAQT LSPFVAVDAL NEALDSYQQV LLTHYGQRM R QKLGFMTEQK EDNALLNELF SLMARERSDY TRTFRMLSLT EQHSAASPLR DEFIDRAAFD DWFARYRRRL QQDEVSDSER QQLMQSVNPA LVLRNWLAQR AIEAAEKGDM TELHRLHEAL RNPFSRDRDD YVSRPPDWGK RLEVSCSS
Source	Baculovirus
Target Names	ydiU
Protein Names	Recommended name: UPF0061 protein ydiU
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.