



Recombinant Escherichia coli O127:H6 Protein syd (syd)

Product Code	CSB-EP486448EOB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7UHK9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MDDLTAQALK DFTARYCDAW HEEHKSWPLS EELYGVSPSPC IISTTEDAVY WQPQPFTGEQ NVNAVERAFD IVIQPTIHTF YTTQFAGDMH AQFGDIKLT LQTWSEDDFR RVQENLIGHL VTQKRLKLP TLFIATLEEE LEVISVCNLS GEVCKETLGT RKRTHLASNL AEFLNQLKPL L
Source	E.coli
Target Names	syd
Protein Names	Recommended name: Protein syd
Expression Region	1-181
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.