



# Recombinant Escherichia coli O127:H6 DNA ligase B (ligB)

<b>Product Code</b>	CSB-MP486491EOB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7UM69
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKVWMAILIS ILCWQSSSWA VCPAWSPARA QEEISRLQQQ IKQWDDDDYWK EGKSEVEDGV YDQLSARLTQ WQRCFGEETH RDAMMPPLNG AVIHPVAHTG VRKMADKIAL SLWMRERSDL WVQPKVDGVA VTLVYRDGKL NKAISRGNGL KGEDWTQKVR LISAVPQTVS GPLANSTLQG EIFLKREGHI QQQMGGINAR AKVADLMMRQ DDSDTLNSLG VFWVAWPDGP QLMTDRLKEL ATAGFTLTQT YTRAVKNADE VARVRNAWWK TKLPFVTDGV VVRAAKEPES RHWLPGQAEW LVAWKYQPVA QVAEVKAIQF AVGKSGKISV VASLAPVMLD DKKVQRVNIG SVRRWQEWNI APGDQILVSL AGQGIPRIDD VVWRGAERTK PTPPENRFNS LTCYFASDVC QEQFISRLVW LGSKQVLGLD GIGEAGWRSL HQTHRFEHIF SWLLLTPEQL QNTPGIAKSK SAQLWHQFNL ARQQPFTRWV MAMGIPLTRA ALNASDERSW SLLLLSTEQF WQQLPGTGSG RVRQVIEWKE NAQIKKLGSW LAAQQITGFE P
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ligB
<b>Protein Names</b>	Recommended name: DNA ligase B EC= 6.5.1.2 Alternative name(s): Polydeoxyribonucleotide synthase [NAD(+)] B
<b>Expression Region</b>	1-561
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.