



# Recombinant Escherichia coli O127:H6 Isocitrate dehydrogenase kinase/phosphatase (aceK)

<b>Product Code</b>	CSB-EP486509EOB-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7UPG7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPRGLELLIA QTILQGFDAQ YGRFLEVTSQ AQRFEQADW HAVQQAMKNR IHLYDHHVGL VVEQLRCITN GQSTDA AFL RVKEHYTRLL PDYPRFEIAE SFFNSVYCRL FDHRSLTPER LFIFSSQPER RFRTIPRPLA KDFHPDHGWE SLLMRVISDL PLRLRWQNKSRDIHYIVRHL TETLGTDNLA ESHLQVANEL FYRNKAAWL V GKLITPSGTL PFLPIHQTD DGELFIDTCL TTTAEASIVF GFARSYFMVY APLPAALVEW LREILPGKTT AELYMAIGCQ KHAKTESYRE YLVYLQGCNE QFIEAPGIRG MVMLVFTLPG FDRVFKVIKD RFAPQKEMSA AHVRACYQLV KEHDRVGRMA DTQEFENFVL EKRHISPALM ALLLQEA AEK ITDLGEQIVI RHLYIERRMV PLNIWLEQVE GQQLRDAIEE YGN AIRQLAA ANIFPGDMLF KNFGVTRHGR VVFYDYDEIC YMTEVNFRDI PLPRYPEDEL ASEPWYSVSP GDVFPEEFRH WLCADPRIGP LFEEMHADLF RADYWRALQN RIREGHVEDV YAYRRRQRFS VRFV
<b>Source</b>	E.coli
<b>Target Names</b>	aceK
<b>Protein Names</b>	Recommended name: Isocitrate dehydrogenase kinase/phosphatase Short name= IDH kinase/phosphatase Short name= IDHK/P EC= 2.7.11.5 EC= 3.1.3.-
<b>Expression Region</b>	1-574
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.