



Recombinant Escherichia coli O127:H6 Periplasmic trehalase (treA)

Product Code	CSB-MP486514EOB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7UQ86
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Purity	>85% (SDS-PAGE)
Sequence	EETSVTPQPP DILLGPLFND VQNAKLFDPD KTFADAVPNS YPLMILADYR MQQNQSGFDL RHFVNVNFTL PKEGEKYVPP EGQSLREHID GLWPVLTRST ENTEKWDSLL PLPKPYVVPD GRFREYVYWD SYFTMLGLAE SGHWDKVADM VANFAHEIDT YGHIPNGNRS YYLSRSQPPF FALMVELLAQ HEGDAALKQY LPQMKEYAY WMDGVENLQA GQQEKRIVKL QDGTLLNRYW DDRDTPRPES WVEDIATAKS NPNRPATEIY RDLRSAASG WDFSSRWMDN PQQNLNLTTRT SIVPVDLNSL MFKMEKILSR ASKAAGDNAM ANQYETLANA RQKGIKYLW NDQQGWYADY DLKSHKVRNQ LTAAALFPLY VNAAAKDRAS KMATATKTHL LQPGGLNTTS VKSGQQWDAP NGWAPLQWVA TEGLQNYGQK EVAMDISWHF LTNVQHTYDR EKKLMEKYDV STTGTGGGGG EYPLQDGFQW TNGVTLKMLD LICPKEQPCD NVPATRPLSE STTQPLKQKE AEPTP
Source	Mammalian cell
Target Names	treA
Protein Names	Recommended name: Periplasmic trehalase EC= 3.2.1.28 Alternative name(s): Alpha,alpha-trehalase Alpha,alpha-trehalose glucohydrolase
Expression Region	31-565
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.