



# Recombinant Escherichia coli O127:H6 UPF0061 protein ydiU (ydiU)

<b>Product Code</b>	CSB-MP486526EOB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7US45
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTLSFITRWR DELPETYTAL SPTPLNNARL IWHNTELANT LSIPSSLFKN GAGVWGGENL LPGMSPLAQV YSGHQFGVWA GQLGDGRGIL LGEQLLADGT TMDWHLKGAG LTPYSRMGDG RAVLRSTIRE SLASEAMHYL GIPTRALSI VTSDSPVYRE TVESGAMLMR VAPSHLRFGH FEHFYRREP EKVRQLADFA IRHYWSHLDD EEDKYRLWFT DVVARTASLI AQWQTVGFAG GVMNTDNMSL LGLTLDYGPFGFLDDYEPGF ICNHSDHQGR YSFDNQPAVA LWNLQRLAQT LSPFVAVDAL NEALDSYQQV LLTHYGQRM R QKLGFMTEQK EDNALLNELF SLMARERSDY TRTFRMLSLT EQHSAASPLR DEFIDRAAFD DWFARYRVRL QQDEVTDSE R QQLMQSVNPA LVLRNWLAQR AIEAAEKGDM TELHRLHEAL RNPFSDRDDD YVSRPPDWGK RLEVSCSS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ydiU
<b>Protein Names</b>	Recommended name: UPF0061 protein ydiU
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.