



Recombinant Escherichia coli O127:H6 Deoxycytidine triphosphate deaminase (dcd)

Product Code	CSB-EP486533EOB-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7UTA3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Purity	>85% (SDS-PAGE)
Sequence	MRLCDRDIEA WLDEGRSLIN PRPPVERING ATVDVRLGNK FRTFRGHTAA FIDLSGPKDE VSAALDRVMS DEIVLDEGEA FYLHPGELAL AVTLESVTLP ADLVGWLDGR SSLARLGLMV HVTahrIDPG WSGCIVLEFY NSGKLPLALR PGMLIGALSF EPLSGPAARP YNRREDAKYR NQQGAVASRI DKD
Source	E.coli
Target Names	dcd
Protein Names	Recommended name: Deoxycytidine triphosphate deaminase Short name= dCTP deaminase EC= 3.5.4.13
Expression Region	1-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.