



Recombinant Escherichia fergusonii Altronate oxidoreductase (uxaB)

Product Code	CSB-BP486745EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LRD0
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MKTLNRRDFP GAQYPERIIQ FGEGNFLRAF VDWQIDLLNE HTDLNSGVVV VRPIETSFPP SLSTQDGLYT TIIRGLNEKG EAVSDARLIR SVNREISVYS EYDEFLKLAH NPEMRFVFSN TTEAGISYHA GDKFDDAPAV SYPAKLTRLL FERFSHFNGA QDKGWIIIPC ELIDYNGDAL RELVLYAQE WALPEAFIQW LDQANSFCST LVDRIVTGYR RDEVTKLEEE LGYHGDGFLDT AEHFYLFVIQ GPKSLATELR LDKYPLNVLV VDDIKPYKER KVAILNGAHT ALVPVAFQAG LDTVGEAMND AEICAFVEKA IYEEIIPVLD LPREELESFA SAVTGRFRNP YIKHQLLSIA LNGMTKFRTR ILPQLLAGQK ANGTLPARLT FALAALIAFY RGERNGETYP VQDDAHWLER YQQLWSQHRD GVIGTQELVA IVLAEKDHWE QDLTQVPLV ELVANDLDAI LEKGMREAVR PLC
Source	Baculovirus
Target Names	uxaB
Protein Names	Recommended name: Altronate oxidoreductase EC= 1.1.1.58 Alternative name(s): Tagaturonate dehydrogenase Tagaturonate reductase
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.