



Recombinant *Escherichia fergusonii* D-amino acid dehydrogenase small subunit (dadA)

Product Code	CSB-EP486756EOR-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LSJ5
Product Type	Recombinant Protein
Immunogen Species	<i>Escherichia fergusonii</i> (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MRVVILGSGV VGVTSAWYLS QAGHDVTVID REPGALETSAANAGQISPG YAAPWAAPGV PLKAIKWMFQ RHAPLAVRLD GTQFQLKMMW QMLRNCDTSH YMENKGRMVR LAEYSRDCLK ELRAATGIEY EGRQGGTLQL FRTEQQFENA TRDIAVLEDA GVPYQLLESN RLAEVEPALA AVAHKLTGGL RLPNDETGDC QLFTQRLAQM AEQAGVKFRF NTPVDKLLFE GEQIYGVKCG DEVIKADAYV MAFGSYSTAM LKGIVDIPVY PLKGYSLTIP VANEDGAPVS TILDETYKIA ITRFDNRIRV GGMAEIVGFN TELLQPRRET LEMVVRDLYP RGGHIEQATF WTGLRPMPDTPD GTPVVGRTSF KNLWLNTGHG TLGWTMACGS GQLLSDILSG RTPAIPYDDL SVARYSSGFS PARPRHLHGV HN
Source	E.coli
Target Names	dadA
Protein Names	Recommended name: D-amino acid dehydrogenase small subunit EC= 1.4.99.1
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.