



Recombinant Escherichia fergusonii Glucans biosynthesis protein G (mdoG)

Product Code	CSB-YP486760EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LT86
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	FSIDDVAK QAQSLAGKGY EAPKSNLPSV FRDMKYADYQ QIQFNSDKAY WSNLKTPFKL EFYHQGMYFD TPKINEVTA TTVKRIKYSY DYFNFGDVKH DKDTVKDLGF AGFKVLYPIN SKDKNDEILS MLGASYFRVI GAGQVYGLSA RGLAIDTALP SGEEFPRFRE FWIERPKPTD KRLTIYALLD SPRATGAYRF VVMPPGRDTVV DVQSKVYLRD KVGKLGVA PL TSMFLFGPNQ PSPATNYRPE LHDSNGLSIH AGNGEWIWRP LNNPKHLAIS SYAMENPQGF GLLQRGRDFS RFEDLDDRYD LRPSAWVTPK GEWGKGVKVEL VEIPTNDET NDNIVAYWTPD QLPEPGKEMN FKYTITFSRD EDKLHAPDNA WVQQTRRSTG DVKQSNLIRQ PDGTIAFVVD FTGTEMKKLP QDTPVTAQTS IGDNGEIVES SVRYNPAIKG WRLTMRVKVK DAKKPTMRA ALATADKTLS ETWSYQLPAN E
Source	Yeast
Target Names	mdoG
Protein Names	Recommended name: Glucans biosynthesis protein G
Expression Region	23-511
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.