



# Recombinant Escherichia coli O8 Ribosomal RNA small subunit methyltransferase H (rsmH)

<b>Product Code</b>	CSB-BP486811E00
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7M125
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O8 (strain IA11)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMENYKHTTV LLDEAVNGLN IRPDGIYIDG TFGRGGHSRL ILSQLGEEGR LLAIDRDPQA IAVAKTIDDP RFSIIHGPF S ALGEYVAERD LIGKIDGILL DLGVSSPQLD DAERGF SFMR DGPLDMRMDP TRGQSAEAWL QTAAEADIAW VLKTYGEERF AKRIARAIVE RNREQPMTRT KELAEEVAAA TPVKDKFKHP ATRTFQAVRI WVNSELEEIE QALKSSLNVL APGGRLSIIS FHSLEDRIK RFMRENSRGP QVPAGLPMTE EQLKKLGGRRQ LRALGKLMPG EEEVAENPRA RSSVLRIAER TNA
<b>Source</b>	Baculovirus
<b>Target Names</b>	rsmH
<b>Protein Names</b>	Recommended name: Ribosomal RNA small subunit methyltransferase H EC=2.1.1.199 Alternative name(s): 16S rRNA m(4)C1402 methyltransferase rRNA (cytosine-N(4)-)-methyltransferase RsmH
<b>Expression Region</b>	1-313
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.