



# Recombinant Escherichia coli O45:K1 Threonine--tRNA ligase (thrS)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP486891EOK   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B7MAS9  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O45:K1 (strain S88 / ExPEC)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MPVITLPDGS QRHYDHAVSP MDVALDIGPG LAKACIAGRV NGELVDACDL<br>IENDAQLSII TAKDEEGLEI IRHSCAHLG HAIKQLWPHT KMAIGPVIDN<br>GFYYDVDLDR TLTQEDVEAL EKRMHELAEK NYDVIKKKVS WHEARETFAN<br>RGESYKVSIL DENIAHDDKP GLYFHEEYVD MCRGPHVPMN RFCHHFKLMK<br>TAGAYWRGDS NNKMLQRIYG TAWADKKALN AYLQRLEEEA KRDRKIGKQ<br>LDLYHMQEEA PGMVFWHNDG WTIFRELEVF VRSKLKEYQY<br>QEVKGPFMMD RVLWEKTGHW DNYKDAMFTT SSENREYCIK<br>PMNCPGHVQI FNQGLKSYRD LPLRMAEFGS CHRNEPSGSL HGLMRVRGFT<br>QDDAHIFCTE EQIRDEVNGC IRLVYDMYST FGFEKIVVKL STRPEKRIGS<br>DEMWDRAEAD LAVALEENNI PFEYQLGEGA FYGPKIEFTL YDCLDRAWQC<br>GTVQLDFSLP SRLSASYVGE DNERKVPVMI HRAILGSMER FIGILTEEFA<br>GFFPTWLAPV QVIMNITDS QSEYVNETQ KLSNAGIRVK ADLRNEKIGF<br>KIREHTLRRV PYMLVCGDKE VESGKVAVRT RRGKDLGSMV VNEVIEKLQQ<br>EIRSRSLKQL EE |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | thrS  |
| <b>Protein Names</b>     | Recommended name: Threonine--tRNA ligase EC= 6.1.1.3 Alternative name(s):<br>Threonyl-tRNA synthetase Short name= ThrRS   |
| <b>Expression Region</b> | 1-642   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.