



# Recombinant Escherichia coli O45:K1 L-rhamnonate dehydratase (rhmD)

<b>Product Code</b>	CSB-MP486927EOK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7MG14
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O45:K1 (strain S88 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MENIMTLPKI KQVRAWFTGG AIAEKGAGGG DYHDQGANHW IDDHIATPMS KYRDYEQSRQ SFGINVLGTL VVEVEAENGQ TGFAVSTAGE MGCFFIVEKHL NRFIEGKCVS DIKLIHDQML NATLYYSGSG GLVMNTISCV DLALWDLFGK VVGLPVYKLL GGAVRDEIQF YATGARPDLA KEMGFIGGKM PTHWGPHDGD AGIRKDAAMV ADMREKCGED FWLMLDCWMS QDVNYATKLA HACAPYNLKW IEECLPPQQY EGYRELKHNA PAGMMVTSGE HHGTLQSFRT LSETGIDIMQ PDVGWCGGLT TLVEIAIAK SRGQLVPHG SSVYSHHAVI TFTNTPFSEF LMTSPDCSTM RPQFDPILLN EPVPVNGRIH KSVLDKPGFG VELNRDCNLK RPYSH
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rhmD
<b>Protein Names</b>	Recommended name: L-rhamnonate dehydratase Short name= RhamD EC= 4.2.1.90
<b>Expression Region</b>	1-405
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.